


GRID CABINET

DESIGNED FOR USE IN RADIOGRAPHIC TABLES AND WALL STANDS

The QUEST Grid Cabinet (Q-1717GC) is designed to fit into most radiographic tables and wallstands that currently house a standard style grid cabinet.

This grid cabinet accepts a variety of X-Ray grids and most ion chambers.



Shown with
QUEST Cassette Tray (Q-1717CT)
& Cassette Sold Separately



Shown with
QUEST Cassette Tray (Q-1717CT)
Sold Separately

KEY FEATURES

CAT# Q-1717GC

GRID CABINET

- Tray locks provide for safe placement and removal of detector
- Accepts X-Ray grids (of 17 1/4" x 18 7/8")
- Accepts ion chambers (of 18" x 18 1/2")
- For Radiographic Tables and Wall Stands

**Specifications are subject to change*



PROVIDES ACCURATE POSITIONING FOR
PORTRAIT & LANDSCAPE ALIGNMENT

WWW.QUESTXRAY.COM

